

IHC: Staining Protocol - Free Floating

This protocol is suitable for the immunohistochemical analysis of formaldehyde or glyoxal fixed vibratome- and cryo-tissue-sections. The tissue-sections are stained free floating. For tissue preparation please refer to our [tissue preparation protocol - formaldehyde](#) or [tissue preparation protocol - glyoxal](#).

Important: For good detection, some proteins have special tissue preparation or [antigen retrieval \(AGR\)](#) requirements. Please refer to the **remarks** section for immunohistochemistry (IHC) on the respective data sheet.

Materials and reagents

- **TBS:** Tris buffered saline, 50 mM Tris pH 7.2/150 mM NaCl
- **Blocking buffer:** 10% normal serum, 0.3% Triton X-100 in TBS (normal serum from the host-species of the secondary antibodies is recommended for blocking).
- **Incubation buffer:** 5% normal serum, 0.3% Triton X-100 in TBS (normal serum from the host-species of the secondary antibodies is recommended for blocking).
- Mounting medium
- **Optional:** DAPI nuclear stain

Procedure

1. Transfer the free floating sections into a staining dish containing **TBS**.
2. **Optional:** Perform [antigen retrieval \(AGR\)](#).
3. Transfer the sections to the **blocking buffer** and block for 1 h at room temperature (RT) (orbital shaker: 70-80 rpm).
4. Transfer the sections to the **incubation buffer** with the primary antibody.
5. Incubate overnight at 4°C (orbital shaker: 60 rpm). **Important:** Some antibodies require incubation at RT. Please refer to the corresponding antibody data sheet.
6. Wash three times for 10 min in **TBS** (RT; orbital shaker: 70 - 80 rpm).
7. Transfer the sections to the **incubation buffer** with the secondary antibody diluted to the manufacturer's recommended concentration.
8. Incubate for 1 h at RT (orbital shaker: 70 - 80 rpm). **Important Notes:**
 - This step can be omitted when fluorophore conjugated primary antibodies are used.
 - In **multiplex** staining, make sure to use secondary antibodies cross-adsorbed against the host species of the other primary antibody used in your experiment. Ideally, all secondary antibodies should come from the same host species. If not, make sure that they have been cross-adsorbed against IgGs of the host-species of the other secondary antibody as well. This avoids cross-reaction between the secondary antibodies.
 - Avoid bright light when working with the secondary antibody to minimize photo bleaching of the fluorescent dye.
9. Wash three times for 10 min in **TBS** (RT; orbital shaker: 70 - 80 rpm). **Optional:** Add DAPI to the first TBS washing step.
10. Wash sections with tap water.
11. Mount slides and microscope.

Note: The SYSY standard protocol generates good staining results in the SYSY labs and may be used as suggestion. However, to achieve the highest specific signal and lowest non-specific background signal, the best antibody concentration, incubation temperature and incubation time must be determined individually.