

Rudolf-Wissell-Str. 28a 37079 Göttingen, Germany

Phone: +49 551-50556-0 +49 551-50556-384 Fax: E-mail: sales@sysy.com Web: www.sysy.com

# VGAT (SLC32A1) cytoplasmic

# domain

Cat.No. 131 003; Polyclonal rabbit antibody, 50 µg specific antibody (lyophilized)

# **Data Sheet**

Reconstitution/ Storage	50 $\mu g$ specific antibody, lyophilized. Affinity purified with the immunogen. Albumin and azide were added for stabilization. For <b>reconstitution</b> add 50 $\mu$ l H <sub>2</sub> O to get a 1mg/ml solution in PBS. Then aliquot and store at -20°C to -80°C until use. Antibodies should be stored at +4°C when still lyophilized. Do not freeze! For detailed information, see back of the data sheet.
Applications	WB: 1: 1000 up to 1: 5000 (AP staining) IP: yes ICC: 1: 100 up to 1: 1000 IHC: 1: 500 up to 1: 3000 IHC-P: 1: 200 up to 1: 500 ExM: external data (see remarks) EM: external data (see remarks)
Immunogen	Synthetic peptide corresponding to residues near the amino terminus of rat VGAT (UniProt Id: O35458)
Reactivity	Reacts with: human (Q9H598), rat (O35458), mouse (O35633), monkey. Other species not tested yet.
Specificity	K.O. validated PubMed: 19052203
Matching control	131-0P
Remarks	<ul> <li>WB: To avoid protein aggregation, do not heat samples for SDS-PAGE.</li> <li>ExM: This antibody has been successfully applied and published for this method by customers (see application-specific references).</li> <li>EM: This antibody has been successfully applied and published for this method by customers (see application-specific references).</li> </ul>

TO BE USED IN VITRO / FOR RESEARCH ONLY NOT TOXIC, NOT HAZARDOUS, NOT INFECTIOUS, NOT CONTAGIOUS

#### Background

The vesicular GABA transporter VGAT is responsible for uptake and storage of GABA and glycine by synaptic vesicles in the central nervous system. For this reason it is frequently referred to as the  $\mathbf{v}$ esicular inhibitory aminoacid transporter VIAAT. It is different from the plasma membrane transporters in that it is driven by a proton electrochemical gradient across the vesicle membrane. So far, only one isoform is known. VGAT is currently the best marker for inhibitory nerve terminals.

#### Selected References for 131 003

Unique luminal localization of VGAT-C terminus allows for selective labeling of active cortical GABAergic synapses. Martens H. Weston MC. Boulland JL. Grønborg M. Grosche J. Kacza J. Hoffmann A. Matteoli M. Takamori S. Harkany T. Chaudhry FA, et al.

The Journal of neuroscience: the official journal of the Society for Neuroscience (2008) 2849: 13125-31.. WB, ICC, IHC; KO verified; tested species: rat

The complement inhibitor CD59 is required for GABAergic synaptic transmission in the dentate gyrus. Wen L. Yang X. Wu Z. Fu S. Zhan Y. Chen Z. Bi D. Shen Y Cell reports (2023) 424: 112349. . WB, IHC, FACS; tested species: mouse

A hypothalamic dopamine locus for psychostimulant-induced hyperlocomotion in mice.

Korchynska S, Rebernik P, Pende M, Boi L, Alpár A, Tasan R, Becker K, Balueva K, Saghafi S, Wulff P, Horvath TL, et al. Nature communications (2022) 131: 5944. . EXM, EM; tested species: mouse

Distribution of gephyrin-immunoreactivity in the trigeminal motor nucleus: an immunohistochemical study in rats. Li Z, Ge S, Zhang F, Zhang T, Mizuno N, Hioki H, Kaneko T, Gao G, Li J

Anatomical record (Hoboken, N.J.: 2007) (2012) 2954: 641-51. . IHC, EM; tested species: rat

Nr4a1 regulates cell-specific transcriptional programs in inhibitory GABAergic interneurons. Huang M, Pieraut S, Cao J, de Souza Polli F, Roncace V, Shen G, Ramos-Medina C, Lee H, Maximov A Neuron (2024):.. IHC, ICC; tested species: mouse

GABAergic neuronal lineage development determines clinically actionable targets in diffuse hemispheric glioma, H3G34-mutant.

Liu I, Alencastro Veiga Cruzeiro G, Bjerke L, Rogers RF, Grabovska Y, Beck A, Mackay A, Barron T, Hack OA, Quezada MA, Molinari V, et al.

Cancer cell (2024):.. ICC, WB; tested species: human

TRPM2 and CaMKII Signaling Drives Excessive GABAergic Synaptic Inhibition Following Ischemia. Burch AM, Garcia JD, O'Leary H, Haas A, Orfila JE, Tiemeier E, Chalmers N, Smith KR, Quillinan N, Herson PS The Journal of neuroscience: the official journal of the Society for Neuroscience (2024) 4419: ... ICC, IHC; tested species: mouse

The TMEM132B-GABAA receptor complex controls alcohol actions in the brain. Wang G. Peng S. Reves Mendez M. Keramidas A. Castellano D. Wu K. Han W. Tian O. Dong L. Li Y. Lu W. et al. Cell (2024) 18723: 6649-6668.e35. . WB, ICC; tested species: mouse

The human milk component myo-inositol promotes neuronal connectivity.

Paquette AF, Carbone BE, Vogel S, Israel E, Maria SD, Patil NP, Sah S, Chowdhury D, Kondratiuk I, Labhart B, Morrow AL, et al. Proceedings of the National Academy of Sciences of the United States of America (2023) 12030: e2221413120. . ICC, IHC; tested species: rat

Enriched environment attenuates hippocampal theta and gamma rhythms dysfunction in chronic cerebral hypoperfusion via improving imbalanced neural afferent levels.

Zheng J. Peng S. Cui L. Liu X. Li T. Zhao Z. Li Y. Hu Y. Zhang M. Xu L. Zhang J. et al. Frontiers in cellular neuroscience (2023) 17: 985246. . WB, IHC; tested species: rat

Spinal Cord Neuronal Network Formation in a 3D Printed Reinforced Matrix-A Model System to Study Disease Mechanisms.

Fischhaber N, Faber J, Bakirci E, Dalton PD, Budday S, Villmann C, Schaefer N

Advanced healthcare materials (2021) 1019: e2100830. . WB, ICC; tested species: mouse

Access the online factsheet including applicable protocols at https://sysv.com/product/131003 or scan the OR-code.



# FAQ - How should I store my antibody?

# **Shipping Conditions**

 All our antibodies and control proteins / peptides are shipped lyophilized (vacuum freezedried) and are stable in this form without loss of quality at ambient temperatures for several weeks.

## Storage of Sealed Vials after Delivery

- Unlabeled and biotin-labeled antibodies and control proteins should be stored at 4°C before reconstitution. They must not be stored in the freezer when still lyophilized!
   Temperatures below zero may cause loss of performance.
- Fluorescence-labeled antibodies should be reconstituted immediately upon receipt. Long term storage (several months) may lead to aggregation.
- **Control peptides** should be kept at -20°C before reconstitution.

# Long Term Storage after Reconstitution (General Considerations)

- The storage freezer must not be of the frost-free variety ("no-frost freezer"). This cycle
  between freezing and thawing (to reduce frost-build-up), which is exactly what should be
  avoided. For the same reason, antibody vials should be placed in an area of the freezer that
  has minimal temperature fluctuations, for instance towards the back rather than on a door
  shelf.
- Aliquot the antibody and store frozen (-20°C to -80°C). Avoid very small aliquots (below 20 µl)
  and use the smallest storage vial or tube possible. The smaller the aliquot, the more the stock
  concentration is affected by evaporation and adsorption of the antibody to the surface of the
  storage vial or tube. Adsorption of the antibody to the surface leads to a substantial loss of
  activity.
- The addition of glycerol to a final concentration of 50% lowers the freezing point of your stock and keeps your antibody at -20°C in liquid state. This efficiently avoids freeze and thaw cycles.

# **Product Specific Hints for Storage**

### Control proteins / peptides

• Store at -20°C to -80°C.

#### **Monoclonal Antibodies**

- Ascites and hybridoma supernatant should be stored at -20°C up to -80°C. Prolonged storage at 4°C is not recommended! Unlike serum, ascites may contain proteases that will degrade the antibodies.
- **Purified IgG** should be stored at -20°C up to -80°C. Adding a carrier protein like BSA will increase long term stability. Many of our antibodies already contain carrier proteins. Please refer to the data-sheet for detailed information.

#### **Polyclonal Antibodies**

- Crude antisera: With anti-microbials added, they may be stored at 4°C. However, frozen storage (-20°C up to -80°C) is preferable.
- Affinity purified antibodies: Less robust than antisera. Storage at -20°C up to -80°C is
  recommended. Adding a carrier protein like BSA will increase long term stability. Most of our
  antibodies already contain carrier proteins. Please refer to the data-sheet for detailed
  information.

#### Fluorescence-labeled Antibodies

• Store as a liquid with 1:1 (v/v) glycerol at -20°C. Protect these antibodies from light exposure.

# Avoid repeated freeze-thaw cycles for all antibodies!

# FAQ - How should I reconstitute my antibody?

#### Reconstitution

- All our purified antibodies are lyophilized from PBS. To reconstitute the antibody in PBS, add
  the amount of deionized water given in the respective datasheet. If higher volumes are
  preferred, add water as mentioned above and then the desired amount of PBS and a
  stabilizing carrier protein (e.g. BSA) to a final concentration of 2%. Some of our antibodies
  already contain albumin. Take this into account when adding more carrier protein.
   For complete reconstitution, carefully remove the lid. After adding water, briefly vortex the
  solution. You can spin down the liquid by placing the vial into a 50 ml centrifugation tube filled
  with paper.
- If desired, add small amounts of azide or thimerosal to prevent microbial growth. This is especially recommended if you want to keep an aliquot a 4°C.
- After reconstitution of fluorescence-labeled antibodies, add 1:1 (v/v) glycerol to a final
  concentration of 50%. This lowers the freezing point of your stock and keeps your antibody in
  liquid state at -20°C.
- Glycerol may also be added to unlabeled primary antibodies. It is a suitable way to avoid freezethaw cycles.
- Please refer to our tips and hints for subsequent storage of reconstituted antibodies and control peptides and proteins.